

Experience the European Design Difference



Operator Furniture

Experience the Difference in Quality

Much More Than Generic Casework

European Design manufactures a wide array of dental operatory furniture to help you and your staff work more efficiently. The European Design Difference is visible in everything we manufacture. Our furniture is a unique blend of concept and construction specifically designed to meet your unique storage, workflow, infection control and esthetic requirements.

Superb Quality for an Excellent Value

Our company was built on the premise of bringing long-term investment value to our customers. We understand that your decision to purchase dental furniture will impact your practice for up to 30 years. If you are looking for superb quality at a great value, then look no further. European Design will exceed your expectations with over 30 years of experience and thousands of dental office installations nationwide. We stand behind our pledge for quality because our furnishings must stand behind you.



Ability to Customize to Meet Your Design Requirements

We offer a variety of sizes and styles to suit your needs. If your practice requires a specialized design or size to meet your specifications, European Design welcomes the opportunity. Just contact us to consult with one of our experienced representatives who will be happy to assist you with your dental furniture requirements.

“Not only are their furnishings excellent looking, they are function-smart. Everything is right where I need it. European Design worked with me to customize the furniture according to my layout needs.”

Christopher A. Troscinski, DDS
Woodridge, Illinois



Exceptional Features Standard in Every Product

You're ready to make a furniture purchase, and there are many choices. When you are comparing features, you can rest-assured that with European Design's dental furnishings, the best features are not optional, they are standard.

The best is not optional, it's standard! Enjoy the advantage.

Features	Local Cabinetmaker	Other Manufacturers	European Design The Right Choice
Seamless Stainless Steel Bases and Support Columns	Not Available Wood/laminate base	Not Available Riveted steel or wood laminate	STANDARD
Solid Nickel Asepsis Handles	Not Available Plastic or aluminum pulls	Not Available Tubular shape pulls (difficult to clean)	STANDARD
Full Extension Ball Bearing Slides	Not Available Partial extension roller-type slides	Not Available Met-a-box 3/4 extension drawer	STANDARD
Soft Close Doors and Drawers	Not Available	Not Available	STANDARD
3/4" Panel Construction	Available (certain situations)	Not Available 5/8" or thinner	STANDARD
3mm Color Matched Edge Banded Doors and Drawers	Not Available Self-edge laminating strips	Available	STANDARD
HPL Laminate or Thermal Fused Melamine Applied Inside and Out	Not Available (Varnished, painted or left unfinished)	Standard	STANDARD
Esthetic Post-Formed Accent Panels	Not Available	Available (from certain companies)	STANDARD
Ability to Customize for Design Requirements	Not Available	Not Available or Discouraged Slows down the manufacturing process	STANDARD

Note - Custom designed cabinets may also be upgraded to include these quality features.

Treatment Stations

Efficient Management of Materials and Workflow

European Design Treatment Stations are designed to enable you to work more efficiently. Materials are managed and meet your ergonomic workflow requirements. Choose between sliding or pivoting type work surfaces and select a built-in handpiece control or nitrous oxide system. Further accessorize your Treatment Station with computer provisions, window inserts, glass shelves or other specialized options.

ET4380 Treatment Station Pivoting L/R Change Work Surface



NEW
Assistant Side Kick with
Pneumatic Open/Close Function



Glass shelves



Equipment drawer

ET4375 Treatment Station
Sliding L/R Change Work Surface



CPU Side Storage

Easy Left/Right Change in Three Quick Steps...



1. Lift small top upward.



2. Slide tops over.



3. Reinsert top on the opposite side.

Central Stations

The Ideal Solution for Multi-room Plans

Communication between you and your patients is confidential and HIPAA guidelines require reasonable steps be taken to ensure this confidentiality. These issues are addressed by providing customized dental furnishings that help to reduce sound and increase privacy in open bay concepts. The Tall Boy Central Station offers you this by extending height and adding upper storage.



RCS104-TB Tall Boy Central Station

The Central Station Tall Boy is the ideal solution for the multi-room plan. Saving you valuable floor space, Tall Boy offers you appropriate storage on either side for both doctor and assistant.

The completely enclosed x-ray compartment can be utilized for two rooms providing you with additional savings.

The Top Trim Panel can be modified to accommodate different ceiling heights.




Esthetic round sink with a solid surface top and an integrated sink



Slide-out doctor's writing surface and drawer storage

In addition to saving floor space, an open bay layout helps you to economize further. This concept eliminates the need for separate doctor and assistant side sink stations and additional x-ray machines.



“During the construction of my new ten operatory dental offices last year, I found the purchase and installation of European Design's furnishings to be the easiest part of getting my practice up and running. Out of all the major brands on the market, European Design was the only one that was able to customize the cabinets to meet my needs. The furnishings are user-friendly, of excellent quality and the service is phenomenal.”

Mark Mueller, DDS
Deer Creek Dental Clinic, Tomah, Wisconsin



Waste chute with easy one hand change bag collar



Concealed glove-tissue-mask dispenser

Round End Side Stations

Providing Better Access to the Operatory

European Design Side Stations feature our unique Round End Storage Compartment or Round Sink End and are available in various designs. Choose to place the sink at the entrance or at the opposite end in reference to your treatment room. Utilize the alcove area of the doctor side station for either a side delivery system or desk area.


Wash stations include a Delta Lever faucet, soap dispenser, waste chute and airless toe-kick control. Temperature and water flow may be directly controlled at faucet.

Side Stations feature soft close doors and drawers that easily and fully close with a touch.



RS5800-R Side Station

The RS5800 Side Station features a curved door and two fixed shelves.




Functional round-end curved door



Solid surface counter tops include an integrated backsplash, a waste chute and choice of integrated or stainless steel undermounted sink.

RS5400-R Side Station

The RS5400 features a round end sink storage compartment.



RS7300-L Doctor Side Station

The RS7300 features a round end shelf storage compartment with a curved door and two fixed shelves.

The alcove area may be used for a side delivery system or desk.

Wall mounted cabinet dispenses gloves, tissues, masks, towels and plastic cups.



The Soft Close Drawers fully extend to reveal the entire contents of drawers. Drawers close smoothly and fully with just a touch.



Stainless steel undermount sink

Specialty Cabinets

Efficient Cabinetry for Your Special Needs

Mobile Cabinets

EML 1820 Mobile

Features a fold-down equipment drawer, left or right side pull-out work surface, drawer storage, floor casters, fixed top and pull handle.



EMF 1820 Mobile

Features side-to-side or front to back sliding top with amalgamator well, drawer storage, all floor casters and pull handle.



Ortho & Pedo



PB-250 Pedo Bench

The Pedo Bench provides a complete pedo work center including drawer storage, 2-handpiece control, HVE, saliva ejector and syringe.



OC1323 Ortho Delivery

Features stainless steel fountain, and panel mount 2-handpiece control system with HVE, saliva ejector and syringe. Note - also available less control system.



OC2932 Ortho Delivery

Features a slide-out work surface, drawer storage, CPU compartment, and panel mount 2-handpiece control system with HVE, saliva ejector and syringe. Optional features: Assortment of monitor and keyboard mounts.

Digital Imaging

SX3236 ScanX® Cabinet

The ScanX® Dental Cabinet is convenient all-in-one storage specifically for use with Scan-X. Features include a slide-out eraser shelf, waste drawer and CPU compartment. Ample space is provided for your LCD display, keyboard and ScanX on the counter surface.




Details

The Details That Add Value to Your Investment



Concealed Hinges -
Lifetime Warranty



Airless Toe-kick Faucet
Control System



Solid Nickel Handles -
Lifetime Warranty



Durable Color Matched
3mm PVC Edgebanding and
Post-formed Accent Panels



Stainless Steel Seamless
Bases with Levelers -
Lifetime Warranty



Full Extension Soft Close
Doors and Drawers -
Lifetime Warranty



3/4" Panel/Dowel and
Dado Construction



U.L. Listed/Hospital
Grade Wiring



Solid-Surface Tops -
10 Year Warranty



Waste Chutes with
Easy-change Bag Collar



European Design Furnishings are designed to withstand rigorous daily use and provide years of trouble free service. Our products utilize some of the finest hardware components available and even go beyond traditional standards offering unique features like our stainless steel bases and solid nickel handles. The integration of high quality dental equipment such as built-in flowmeters and control systems further enhance our operatory furnishings, making it easy to select a complete cabinet system that meets all your requirements.

European Design's network of full service dealers are available to assist you with your practice planning. To locate a dealer near you, please contact us toll-free at 1-800-257-7407.

European Design, Inc.
115 G.L. Comer Road, Glasgow Kentucky 42141
europeandesigninc.com
1-800-257-7407
Fax 270-651-1732

© 2008 Midmark Corporation
60 Vista Drive, P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0696-00 Rev. B1 (9/08)



Because we care.