



PROCENTER WALL & CABINET  
MOUNTED DELIVERY UNITS

UNCOMPROMISING QUALITY

INNOVATIVE DESIGN

## PROCENTER CABINET MOUNTED DELIVERY UNIT

PROCENTER CABINET MOUNTED DELIVERY UNIT, SHOWN WITH THE SHORT MOUNT ARM, 16" FLEX ARM AND OPTIONAL MIDMARK CHAIR CONTROL TOUCHPAD.

### THE MIDMARK PROCENTER

CABINET AND WALL MOUNTED INSTRUMENT DELIVERY SYSTEMS PROVIDE ALL THE ADVANTAGES OF THE CHAIR MOUNTED PROCENTER INCLUDING INDEPENDENT HANDPIECE CIRCUITS FOR OUTSTANDING RELIABILITY, INDIVIDUAL AIR AND WATER NEEDLE VALVE ADJUSTMENTS FOR PRECISE CONTROL, A SELF-CONTAINED WATER SYSTEM TO MANAGE WATER QUALITY AND FACILITATE WATER SYSTEM FLUSHING REGIMENS, AND A CITY WATER / BOTTLE WATER SELECTOR SWITCH. YOUR CHOICE OF 16" OR 24" AIR LOCK FLEX ARMS AND SHORT OR LONG MOUNT ARMS ALLOW YOU TO CUSTOMIZE THESE PROCENTER DELIVERY UNITS TO MEET YOUR OPERATORY DESIGN REQUIREMENTS.

THE LOCATION OF YOUR INSTRUMENT DELIVERY UNIT IS A KEY DECISION IN OPERATORY DESIGN. SIDE DELIVERY OFFERS THE ADVANTAGES OF PLACING THE HANDPIECES IN A COMFORTABLE WORKING POSITION WHILE PROVIDING A FLEXIBLE WORK ZONE. BY MOUNTING THE SIDE DELIVERY UNIT TO A CABINET OR WALL, YOUR INSTRUMENT DELIVERY CAN BE STORED OUT OF VIEW AND ALLOWS FOR UNRESTRICTED CHAIR ACCESS FOR PATIENT ENTRY AND EXIT.



MIDMARK IS AN ISO 9001 CERTIFIED COMPANY

U.S. PATENT No. D458379 AND 6,692,253



FOR MORE INFORMATION OR A DEMONSTRATION, CONTACT YOUR MIDMARK DEALER OR MIDMARK CORPORATION.

MIDMARK CORPORATION,  
60 VISTA DR., P.O. BOX 286  
VERSAILLES, OHIO 45380 U.S.A.  
1-800-MIDMARK  
FAX 1-877-725-6495

OUTSIDE THE U.S.A.  
CALL: 1-937-526-3662  
FAX: 1-937-526-8214

[midmark.com](http://midmark.com)

PROCENTER WALL MOUNTED DELIVERY UNIT, SHOWN WITH THE LONG MOUNT ARM AND 24" FLEX ARM.

## PROCENTER WALL MOUNTED DELIVERY UNIT

© Midmark Corporation 2002  
All products are subject to improvement changes without notice.  
Litho in U.S.A.  
007-0340-00 Rev. A1 (3/06)