

MIDMARK\*



MIDMARK ULTRATRIM TM & ULTRACOMFORT TM CHAIRS

UNCOMPROMISING OUALITY

INNOVATIVE DESIGN

## ULTRATRIM™ DENTAL CHAIR

MIDMARK OFFERS YOU THE CHOICE OF TWO DISTINCTIVE PATIENT CHAIRS, THE ULTRATRIM AND THE ULTRACOMFORT. WHILE VERY DIFFERENT, BOTH CHAIRS SHARE THE SAME MIDMARK DESIGN PRINCIPLES FOR QUALITY DENTAL CHAIRS: COMFORT, ACCESS AND CONTROL.

THE MIDMARK PATENTED CHAIR-BACK DESIGN FEATURES A CAST ALUMINUM BACK WITH INTEGRATED ARMRESTS. THE INTEGRATED ARMRESTS CRADLE THE KEY PATIENT SUPPORT AREAS OF THE ELBOWS AND INNER SHOULDERS DELIVERING A COMFORTABLE, SECURE PATIENT. AND, AS THE CHAIR IS RETURNED TO THE UPRIGHT POSITION, THE ARMRESTS ARE AUTOMATICALLY POSITIONED FOR UNENCUMBERED EASE OF ENTRY AND EXIT.

THE CHAIR-BACK PIVOT AND SYNCHRONIZED MOVEMENT FURTHER ENHANCE THE PATIENT'S COMFORT AND SUPPORT. THE MIDMARK CHAIR-BACK MOVEMENT MATCHES THAT OF THE PATIENT AS THEY MOVE FROM A SLIGHT RECLINED TO SUPINE POSITION, ELIMINATING THE NEED FOR PATIENT OR HEADREST READJUSTMENT.

THE MIDMARK ULTRATRIM DENTAL CHAIR IS THE ULTIMATE IN CLINICIAN ACCESS. WE CAREFULLY PLOTTED ESSENTIAL SUPPORT AREAS FOR A WIDE RANGE OF BODY TYPES, WHICH ALLOWED US TO REMOVE MATERIALS IN JUST THE RIGHT PLACES. THE RESULT? THE ULTRATRIM CHAIRBACK IS OVER 40% THINNER AND NARROWER THAN OUR STANDARD CHAIR-BACK.









THE MIDMARK ULTRACOMFORT PATIENT CHAIR IS THE ULTIMATE IN PATIENT COMFORT. WITH ITS FULL CONTOURED BACK AND SPECIALLY DESIGNED LUMBAR SUPPORT, THE ULTRACOMFORT CHAIR WILL KEEP PATIENTS SUPPORTED AND COMFORTABLE THROUGH EVEN THE LENGTHIEST OF PROCEDURES.

CHAIR-BACK DESIGN IS KEY TO PROPER CLINICIAN POSITIONING AND PATIENT ACCESS. THE ULTRACOMFORT CHAIR-BACK IS THIN AND NARROW, ALLOWING YOU TO WORK IN AN ERGONOMICALLY CORRECT POSITION. THE ULTRATRIM CHAIR-BACK TAKES THIS TO A NEW LEVEL WITH A BACK SO THIN AND SO NARROW THAT EVEN DIRECT PATIENT ACCESS IS EASILY ACHIEVED.

MULTIPLE CHAIR CONTROLS ARE STANDARD, FEATURING FOUR PROGRAMMABLE PRESETS THAT MAKE PLACING THE CHAIR AND PATIENT IN YOUR PREFERRED POSITION EASY... WITH THE TOUCH OF A BUTTON. AND THE HYDRAULIC DRIVE SYSTEM DELIVERS SMOOTH, QUIET, PRECISE OPERATION.

Choice is good, and with the Midmark UltraTrim and UltraComfort dental chairs, we expand your options. Both chairs are available in seamless Silhouette or plush Ultraleather<sup>tm</sup> upholstery, with your choice of double articulating or magnetic headrest.

ULTRACOMFORT™ DENTAL CHAIR

MIDMARK IS AN ISO 9001 CERTIFIED COMPANY.



Because we care.

FOR MORE INFORMATION OR A DEMONSTRATION, CONTACT YOUR MIDMARK DEALER OR MIDMARK CORPORATION.

MIDMARK CORPORATION 60 VISTA DR., P.O. BOX 286 VERSAILLES, OHIO 45380 U.S.A. 1-800-MIDMARK FAX 1-877-725-6495

Outside the U.S.A. Call: 1-937-526-3662 Fax: 1-937-526-8214

midmark.com

© Midmark Corporation 2004 All products are subject to improvement changes without notice. Litho in U.S.A. 007-0543-00 Rev. A1 (1/05)