

S E A T I N G



## Ergonomically Refined

The Midmark Operator Stools are not only engineered to deliver premium styling and exceptional comfort, they are also ergonomically refined to address musculoskeletal pain common in dentistry. Promoting a more balanced posture, the seating system offers a contoured seat, with raised front, to improve weight distribution and overall comfort. Comfort is also enhanced with waterfall front edges that reduce pressure on the back of the thighs. Multiple seat sizes and seat heights, paired with unsurpassed adjustability, allow you to better personalize your seating solution.



## Midmark Dentist's Stool

Seat size: Standard 17.5" x 16.5" or Small 16.25" x 16"

Up to 12 degrees of forward tilt helps to maintain natural S-curve of the back

One lever adjustment for locking seat and backrest tilt

Ratcheted backrest for intuitive height adjustment without knobs

Cylinder height: Standard 16.5" - 22" or Short 15.25" - 18.3"

Upholstery available in Silhouette and optional Ultraleather™

## Midmark Assistant's Stool

Seat size: 16.25" x 16"

Optional ratcheted backrest for intuitive height adjustment without knobs

Ratchet release body support with height adjustment

Integrated footrest, automatically adjusts with seat height

Cylinder height: Standard 21.7" - 27.5" or Short 17.8" - 20.8"

Upholstery available in Silhouette and optional Ultraleather™



Maintaining an ergonomically correct working posture is critical to perpetuating a pain-free career. Midmark operator stools offer a contoured seat with a raised front that naturally positions the body against the backrest, thus promoting a more balanced posture. Tilting the seat forward allows you to open the hip angle over 90 degrees, which helps to maintain the natural lower back curve and reduce disc pressure. The contoured backrest comfortably supports the lower back while the adjustable armrests help to support the elbows, all in an effort to reduce fatigue in the lower back, neck and shoulders.



## Dentist's Armrests

Optional armrests for the dentist's stool offer support for elbows in the neutral shoulder position, in an effort to reduce fatigue of the neck and shoulder muscles. Arm pad adjustments include height, swivel and slide, optimizing positioning of the armrests to better accommodate your personal needs.

**Midmark is an ISO 13485 and ISO 9001 Certified Company.**  
For more information or a demonstration, contact your Midmark  
dealer or Call: 1-800-MIDMARK Fax: 1-877-725-6495  
Outside the U.S.A. Call: 1-937-526-3662 Fax: 1-937-526-8214  
or visit our website at [midmark.com](http://midmark.com)

© 2008 Midmark Corporation  
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286  
Products subject to improvement changes without notice  
Litho in U.S.A. 007-0715-00 Rev. A1 (4/08)



Because we care.