

L E D L I G H T



Seeing is believing.

The Midmark LED Operatory light is a market-driven solution utilizing cutting edge LED technology to provide the ultimate in dental lighting. It provides the true, reduced shadow white light you need to properly match shades, identify details and diagnose tissue. The Midmark LED offers a feathered edge to minimize eye strain caused by sharp edges with high contrast differences. Time-tested, user-friendly features simplify use of the light as well as minimizing fatigue, time and motion. The Midmark LED accomplishes all of this while using only 22 watts of electricity, less than half of the leading dental LED light.



True composite safe setting: Midmark uses filters to remove the 470nm blue light that cures composites.

Lighting Specifications

Operate Mode Color Temperature:	5700 Kelvin
Operate Mode Intensity:	1400, 1950, 2500 Foot Candles 15000, 21000, 27000 Lux
Composite Mode Color Temperature:	4000 Kelvin
Composite Mode Intensity:	2100 Foot Candles; 22000 Lux
Pattern Size:	3.3" x 7.5"

Retrofit Kit

Bring an existing operatory up to date with our easy to install retrofit kits.

For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214
or visit our website at midmark.com

© 2011 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0936-00 Rev. A1 (3/12)

