

ECO-VAC

DRY VACUUM SYSTEMS

"One Green Machine"

**Enables you to care for your patients,
the environment AND your bottom line.**

- **SAVES OVER 50% ON ELECTRICITY COSTS WHEN COMPARED TO SIMILAR VACUUM PUMPS**
- **ZERO WATER USAGE SAVES WATER AND SEWAGE COSTS**
- **DRY TURBINE SECTION USES NO OIL LIKE SOME OTHERS**
- **ELIMINATES MAINTENANCE COSTS AND KEEPS THE OIL VAPORS FROM THE EQUIPMENT ROOM AND ENVIRONMENT**
- **INTEGRATED AMALGAM SEPARATOR TO REMOVE MERCURY WASTE (OPTIONAL)**
- **VERY QUIET OPERATION AT ONLY 62 DECIBELS**

We utilize a Variable Frequency Drive (VFD) to control the vacuum levels. With no suction being used, the motor is running one third the speed and using very little power. Every time a suction hose is being used the motor automatically speeds up to maintain constant suction.

Maintenance free unit with no traps to clean, no oil to check or change. Truly a ZERO maintenance machine.



***TW* TECH WEST INC.**

Manufacturers of Dental Vacuum
and Air Systems

2625 N. Argyle Ave. • Fresno, CA 93727
(559) 291-1650 • (800) 428-7139 • FAX (559) 348-9677

Our EcoVac system provides the advantages of proven technology, durable construction, quiet operation and oil-free performance with the latest in water- and energy-saving features. This is not only important to you as an environmentally responsible business owner, it also prepares your practice for the upcoming change to metered water usage in many areas of the country.



ECO VAC DRY VAC SYSTEMS

FEATURES

- Green design
- Whisper quiet motor
- Energy-efficient operation
- Oil free performance
- Nickel Coated turbine section
- Delivered ready to install
- 5-year/5000-hour warranty

BENEFITS

- Saves enough water each month to fill two swimming pools
- Reduces stress on patients and staff
- Lowers your power bill and minimizes your carbon footprint
- Provides clean dry operation
- Prevents corrosion and moisture damage inside the turbine.
- Requires no new piping to switch from wet to dry vacuum system, potentially saving thousands in installation costs
- Protects your investment and keeps your practice up and running

ECO DRY-VAC SYSTEMS

MODEL NUMBER	MAX USERS	# OF HEADS	TOTAL H.P.	CFM @ 8 HG	TANK GALLONS	VOLTAGE	TOTAL AMPS	BREAKER SIZE	DB(A)	EQUIP. WEIGHT	EQUIP. HEIGHT	EQUIP. WIDTH	EQUIP. DEPTH	SHIPPING WEIGHT	SHIPPING HEIGHT	SHIPPING WIDTH	SHIPPING DEPTH
VPD5S2	5-7	1	2-3	45	14	208/230	7	20	62	172	52	23	20	226	43	43	30
VPD10D2	10-14	2	4-6	90	20	208/230	14	20	62	365	39	42	20	435	43	43	30

* 2 BOXES

Revised 9/09

FIVE-YEAR, 5000 HR. WARRANTY

Exceeds NFPA 99C Guidelines