

THE ROCKY SERIES

- VERY QUIET MOTOR DESIGN ONLY 66 DB-A
- ALL AIR RECEIVER TANKS HAVE INTERNAL EPOXY COATING APPLIED
- THE AIR IS DRYED BEFORE THE TANK TO PREVENT BACTERIA BUILD UP IN THE TANK
- ALL MODELS HAVE DESICCANT AIR DRYING SYSTEMS
- ROCKING PISTON TECHNOLOGY DECREASES NOISE AND INCREASES MOTOR LIFE
- COMPRESSORS ARE AVAILABLE IN DUAL, TRIPLE OR QUAD HEAD CONFIGURATIONS
- TEFLON COATED PISTONS AND SLEEVES INSURE LONG LIFE OF THE COMPRESSOR HEADS
- REMOTE AIR INTAKE FILTERS TRAP DIRT, DUST AND AIR BORN POLLUTANTS
- ENERGY EFFICIENT THERMALLY PROTECTED SPLIT PHASE MOTORS
- MOISTURE INDICATORS STANDARD ON ALL MODELS



***TW* TECH WEST INC.**

Manufacturers of Dental Vacuum
and Air Systems

2625 N. Argyle Ave. • Fresno, CA 93727
(559) 291-1650 • (800) 428-7139 • FAX (559) 348-9677

The Tech West Rocky Series Oilless Air Compressors provide 100% oil free air to the dental office. With our high performance air drying system, we guarantee the clean dry air you expect. The design of the Tech West Rocky Series Air Compressor provides smooth, quiet, vibration free operation at only 65 dB-a noise level. We offer them in twin, triple or quad head configurations to meet most any need of today's Dental Practice.

Tech West will provide full plumbing and electrical design upon request.

THE ROCKY SERIES

FEATURES

Desiccant Air Dryer with .01 micron filtration

Rocking piston design

USA built electrical motors

Integrated carrying handles

Moisture indicators on the air line

BENEFITS

Assures clean, dry air with a dew point of -25 degrees

Reduces noise and increases piston seal life

Gives the motor a voltage range of 200-253 so a buck boost transformer is not needed

Makes for easy moving and installation

Gives you an early indication of moist air prior to moisture getting into your air lines.

THE ROCKY SERIES

MODEL NUMBER	MAX USERS	AIR DRYER	QUIET COVER	# OF HEADS	TOTAL H.P.	CFM	TANK GALLONS	VOLTAGE	TOTAL AMPS	BREAKER SIZE	DB(A)	EQUIP. WEIGHT	EQUIP. HEIGHT	EQUIP. WIDTH	EQUIP. DEPTH	SHIPPING WEIGHT	SHIPPING HEIGHT	SHIPPING WIDTH	SHIPPING DEPTH
ACOR2D1	2-3	YES	NO	2	1.5	5.2	10	115	11	20	68	170	26	26	18	190	37	35	30
ACOR2D1Q	2-3	YES	YES	2	1.5	5.2	10	115	11	20	63	210	31	26	23	230	42	42	31
ACOR2D2	2-3	YES	NO	2	1.5	5.2	10	208/230	8	20	68	170	26	26	18	190	37	35	30
ACOR2D2Q	2-3	YES	YES	2	1.5	5.2	10	208/230	8	20	63	210	31	26	23	230	42	42	31
ACOR3T2	3-5	YES	NO	3	2.25	7.8	20	208/230	12	20	74	265	30	36	23	274	43	42	31
ACOR3T2Q	3-5	YES	YES	3	2.25	7.8	20	208/230	12	20	71	305	35	37	23	314	43	42	31
ACOR4Q2	4-6	YES	NO	4	3	10.4	30	208/230	16	20	75	290	32	41	24	320	43	42	31

Revised 1/08

**ALL OILLESS COMPRESSORS CARRY A FULL THREE YEAR WARRANTY ON OILLESS HEADS
AND A FULL TWO YEAR WARRANTY ON ALL OTHER COMPONENTS**

ti TECH WEST INC.
Manufacturers of Dental Vacuum
and Air Systems

Exceeds NFPA 99C Guidelines



Complies to UL Standards
ANSI/UL 544.

2625 N. Argyle Ave. • Fresno, CA 93727
(559) 291-1650 • (800) 428-7139 • FAX (559) 348-9677

FDA
RECOGNIZED