

T R A C K L I G H T
M O N I T O R







Technology Brought Into View

The Track Light Monitor from Midmark combines the best in dental lighting with unparalleled monitor positioning, allowing you to enhance the patient experience while improving efficiency and case acceptance.

Entertain

Patients relax while enjoying their favorite DVD or TV program. Acting as a distractant, anxiety is reduced without the use of traditional sedation methods.

Inform

Increase case acceptance by sharing digital x-rays, intraoral camera images and patient education software.

Increase referrals

When patients are pleased with an experience, they tell their friends. This allows you to expand your customer base and attract the "ideal patient."

Increase efficiency

When relaxed, the patient tends to move less during long procedures. Reduced patient movement allows you to complete procedures more quickly, improving your productivity.



Chair Side

Ideally positioned for chair side consultation, the monitor can be placed at eye level to the patient for effortless viewing of digital x-rays, intraoral camera images, patient education and entertainment.



Overhead

The counterbalanced flex arm and independent trolleys allow the monitor to be fluidly positioned over the patient while in a supine position. Interference with the operatory light is minimized so the patient and doctor enjoy unobstructed views.



Connectivity

Multimedia connections include VGA, 3.5mm PC audio, Coaxial cable, HDMI, S-Video and RCA composite.

SEAMLESS INTEGRATION

The monitor system conveniently consolidates inputs from PC, TV and DVD. Clutter is reduced and future expandability is simplified.

17" dental grade monitor with remote control

Sealed display, buttons and speakers withstand the use of disinfectants

Colors: 16.7 million

Resolution: SXGA 1280 x 1024

Digital TV tuner with CC: NTSC, ATSC & QAM

HD support: 480i, 480p, 720p, 1080i



Monitor Mount

The optional 18.5" monitor mount will accommodate customer-supplied monitors. Monitors must be 12V, VESA compliant and weigh from 5 – 13 lbs.

EXCEPTIONAL LIGHTING

Where many lights are an afterthought of old technology, the Midmark light pushes lighting quality to the edge. It features an integrated PC board that monitors and adjusts the voltage to ensure color corrected brightness, LED intensity readouts and touchpad controls. All of this helps to deliver the true, reduced shadow white light you need to properly match shades, identify details and diagnose tissue. User-friendly features simplify using the light while minimizing fatigue, time and motion.

The touchpad control offers three intensity adjustments as well as a composite safe setting to prevent premature curing. The automatic motion controlled on/off function can be set on the touchpad with the simple touch of a button.



Lighting Specifications

Operate Mode Color Temperature:	4500, 4750, 5000 Kelvin
Operate Mode Intensity:	1400, 1950, 2500 Foot Candles 15000, 21000, 27000 Lux
Composite Mode Color Temperature:	4000 Kelvin
Composite Mode Intensity:	550 Foot Candles 5900 Lux
Color Rendering Index:	Greater than 90%
Pattern Size:	3.3" x 7.5"



Precise Positioning

The ergonomic design offers complete control with swivel, tilt and rotate, as well as easy access to the touchpad. Tension adjustments enable you to set the tension and motion desired for your specific needs.



Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark
dealer or Call: 1-800-MIDMARK Fax: 1-877-725-6495
Outside the U.S.A. Call: 1-937-526-3662 Fax: 1-937-526-8214
or visit our website at midmark.com

© 2007 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0646-00 Rev. D1 (5/09)



Because we care.